

PERRY BAROMEDICAL CORPORATION
"Excellence In Hyperbaric Medical Facilities"
FAX: 561/840-0398

May 20, 2005

National Baromedical Services
FAX # 803-434-4354

Re: Review of CAVICIDE Disinfectant Cleaner MSDS Sheets

Perry Baromedical must advise against use of CAVICIDE as a cleaner/disinfectant for the PBC hyperbaric chamber systems. The acrylic tube is susceptible to any form of an alcohol-based solution, and per the MSDS sheets, this material contains 400ppm isopropanol alcohol as an active ingredient. As such, it must not be used.

If you have any further questions, please feel free to call.

Regards,

Keith Evans
PBC Engineering

CC: T. Buck
D. Vermillion